Instruction for use

HEINE Finoff

DEUTSCH | ENGLISH | FRANÇAIS | ESPAÑOL | ITALIANO

Produktübersicht
Product overview
Vue d'ensemble du produit
Partes del producto
Informazioni generali sul prodotto

HEINE Finoff

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HEINE Finoff

Please read and follow these instructions for use and keep them for future reference.

Intended Use

The hand-held, internally powered Finoff transilluminator is intended for general, transient and transillumination of intact tissue and, more specifically, for scleral transillumination. The Finoff transilluminator can also be combined with the blue filter for illuminating the eye. The device may only be used by qualified medical personnel and in professional healthcare facilities.

⚠️ For U.S. only: Federal law restricts this device to sale by or on the order of a Physician or Practitioner.

Warnings and Safety Information

⚠️ CAUTION! This symbol indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

📌 NOTE! This symbol indicates valuable advice. This notes are important, but no related to hazardous situations.

Product overview (see fig. ①)

- 1 HEINE Rechargeable handle
- 2 Finoff Transilluminator
- 3 LED illuminant
- 4 Aperture with Cobalt blue filter
- 5 Aperture

Setting up

To set up the Finoff transilluminator, screw the instrument head into the HEINE battery handle or plug it into the HEINE rechargeable handle.

Verify that the lamp voltage complies with the supply voltage of the handle. The coloured marking on the bottom of the lamp shows you the lamp voltage:

White ring = HEINE 2.5 V bulb

only for use with the HEINE BETA battery handle

Black ring = HEINE 3.5 V bulb

only for use with the HEINE BETA rechargeable handle and HEINE wall transformer.

The setting up and operation of the HEINE power supplies are described in a separate instruction of use.

Operation

📌 We recommend that the examination with the HEINE Finoff transilluminator be carried out in a darkened room. The surface of the eye must be anaesthetised before scleral transillumination.

Set the light intensity as low as possible during the examination.

Please note that a high light intensity can be uncomfortable for the patient and may make it difficult to recognise areas of thin tumour growth.

⚠️ The Finoff transilluminator is intended for transient examination at a maximum intensity of < 60 sec. with a 15 min. break until the next application.

Only use the Finoff transilluminator for scleral transillumination with an aperture (4, 5).

To illuminate the eye, fit the aperture with the blue filter (4) and move the instrument as close as necessary to the eye without touching it.

For scleral transillumination, rest the distal end with the aperture (5) flat on the bulbous without exerting pressure.

For transillumination of skin tissue, place the light-emitting surface on the skin.

If you notice a significant drop in brightness immediately after switching on the device, you should charge the rechargeable handle.

Light exposure hazard

⚠️ Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures. The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

The intensity of the light directed into the patient's eye must be reduced to a minimum level. Infants, aphakes and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

⚠️ The duration of exposure during an eye examination must not exceed 28 sec. with local illumination and a distance of 5 mm at maximum intensity.

The duration of exposure with direct, local and unchanged skin contact and maximum intensity must not exceed 6 seconds.

Hygienic reprocessing

The instruction is available:

- at www.heine.com

- in a paper version which you can request from the address listed

Maintenance

Changing the light source

⚠️ Verify that the lamp voltage complies with the supply voltage of the handle. Allow the device to cool down before changing the bulb.

Detach the Finoff from the handle and pull the bulb out of the guide tube of the connector plug holding it at the narrow collar. Insert the new lamp all the way into the lamp guide.

Service

The device has no components serviceable by the end-user.

General Notes

📌 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit www.heine.com.

The expected life cycle, when the device is normal used and the warning and safety information as well as the maintenance instructions are observed, is up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

Note to the user and | or patient:

All serious incidents that occur in connection with the product must be reported to HEINE Optotechnik GmbH & Co. KG and the member state's competent authority.

General Warnings

⚠️ Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

Do not use the device in presence of flammable gases | liquids, or in an oxygen rich environment.

Do not use the device in fire- or explosive risk area.

This product is not allowed to enter or be used in areas with strong magnetic fields e.g. MRI scanners.

Danger due to loose, falling parts. Check that all parts are firmly attached before each use.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources. Repairs shall only be carried out by qualified persons.

Disposal

♻️ The product must be recycled as separated electrical and electronic devices.

Please observe the relevant state-specific disposal regulations.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Spare parts
- Explanation of the used symbols

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Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

	Inside hospitals except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.
Statement for the operational environments	The supply voltage quality should be that of a typical hospital environment. <p>Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.</p>
Performance features of the ME system that have been determined to be essential to the performance	
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle	None

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Finoff, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

EMC compatibility is only ensured if original HEINE spare parts, accessories and power sources are used. The EMC compatibility when using power sources from other manufacturers must be evaluated by the user.

Test	Compliance test levels
RF emissions CISPR11	Group 1 Class B
RF emissions (EN 55011 CISPR 11)	Passed
Conducted emissions (EN 55011 CISPR 11)	Not applicable, because of internally powered mode
Radiated emissions (EN 55011 CISPR 11)	Passed
Harmonic current emissions (IEC 61000-3-2)	Not applicable, because of internally powered mode
Voltage changes, voltage fluctuations and flicker (IEC 61000-3-3)	Not applicable, because of internally powered mode
Immunity	See attached immunity test levels

Immunity test levels		
Test	IEC 60601-1-2 test levels	Compliance test levels
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: ±8 kV <p>Air Discharge: ±2 kV, ±4 kV, ±8 kV, ±15 kV</p>	
Radiated RF EM fields (IEC 61000-4-3)	3 V/m <p>80–2700 MHz <p>80 % AM at 1 kHz</p></p>	
Electrical fast transients bursts (IEC 61000-4-4)	± 2 kV <p>100 kHz repetition frequency</p>	Not applicable, because of internally powered mode
Surges (IEC 61000-4-5)	± 0.5 kV, ± 1 kV	Not applicable, because of internally powered mode
Conducted disturbances induced by RF fields (IEC 61000-4-6)	3 V <p>0.15 MHz – 80 MHz <p>6 V in ISM bands between 0.15 MHz and 80 MHz <p>80 % AM at 1 kHz</p></p></p>	
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse Modulation: 18 Hz; 27 V/m <p>450 MHz, FM: ±5 Hz deviation; 1 kHz sine; 28 V/m <p>710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m <p>810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m <p>1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m <p>2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; <p>5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m</p></p></p></p></p></p>	
Power frequency magnetic fields (IEC 61000-4-8)	30 A/m; 60 Hz	
Voltage dips (IEC 61000-4-11)	0 % U _r ; 0.5 cycle; at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° <p>0 % U_r; 1 cycle and 70 % U_r; 25/30 cycles Single phase: at 0°</p>	Not applicable, because of internally powered mode
Short interruptions (IEC 61000-4-11)	0 % U _r ; 250/300 cycle	

Technical specification	
Environmental conditions for operation	+10 °C to +35 °C <p>30 % to 75 % rel. humidity <p>700 hPa to 1060 hPa</p></p>
Environmental conditions for storage	+5 °C to +45 °C <p>45 % to 80 % rel. humidity <p>500 hPa to 1060 hPa</p></p>
Environmental conditions for transport	-20 °C to +50 °C <p>45 % to 80 % rel. humidity <p>500 hPa to 1060 hPa</p></p>

Spare parts

The corresponding spare parts are listed under the following link: https://www.heine.com/Finoff-transilluminator-accessories-IFU

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