**Instruction for use**  
**HEINE OMEGA 600 |**  
**HEINE OMEGA 600 wired**

# **Instruction for use**

## **HEINE OMEGA 600 |**

## **HEINE OMEGA 600 wired**

<b>DEUTSCH</b>	<b>3</b>
<b>ENGLISH</b>	<b>7</b>
<b>FRANÇAIS</b>	<b>11</b>
<b>ESPAÑOL</b>	<b>15</b>
<b>ITALIANO</b>	<b>19</b>
<b>SVENSKA</b>	<b>23</b>
<b>NEDERLANDS</b>	<b>27</b>
<b>DANSK</b>	<b>31</b>
<b>NORSK</b>	<b>35</b>
<b>SUOMI</b>	<b>39</b>
<b>PORTUGUÊS</b>	<b>43</b>

# HEINE OMEGA 600

## HEINE OMEGA 600 wired

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# HEINE OMEGA 600

## HEINE OMEGA 600 wired

 Please read and follow these instructions for use and keep them for future reference.

### Intended use

The HEINE OMEGA 600 and HEINE OMEGA 600 wired binocular indirect ophthalmoscope is a head worn AC-powered or battery powered device for transient use, containing illumination and viewing optics intended to examine posterior segments of the eyes.

The product must only be used by qualified medical professionals and in medical healthcare facilities.

### For U.S. only:

 **Federal law restricts this device to sale by or on the order of a Physician or Practitioner.**

### Warnings and safety information


 **CAUTION!** This symbol indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** This symbol indicates valuable advice. Notes are important, but not related to hazardous situations.

### Product overview


OMEGA 600

Fig. 1


OMEGA 600 wired

Fig. 2

Optics unit (13) for OMEGA 600 and OMEGA 600 wired

Fig. 3

Rechargeable battery CB1

Fig. 4

- 1 Width adjustment
- 2 Height adjustment
- 3 Padding
- 4 Socket for brightness control
- 5 Adjustment lever
- 6 Brightness control
- 7 Charge status indicator
- 8 Battery compartment
- 9 USB-C socket
- 10 Rechargeable battery CB1
- 11 Power indicator
- 12 Power supply connection
- 13 Optics unit
- 14 Aperture selection lever
- 15 Illumination height adjustment
- 16 Dust cover
- 17 Stereoscopic adjustment lever
- 18 Filter selection lever
- 19 Eyepieces

### Setting up

To put the OMEGA 600 into operation, insert the CB1 battery (10) into the instrument's battery compartment (8) by quickly pressing it in so that it audibly clicks into place. We recommend to charge the battery CB1 (10) of the device fully before first use.

To put the OMEGA 600 wired into operation, connect the power supply unit to a mains plug using the appropriate country-specific primary adapter. Fig. 5 and Fig. 6 describe the approach for changing the plug.


Fig. 5


Fig. 6

 Connect the power supply unit to the mains after checking the mains voltage indicated on the identification plate.

Mount the power supply in such a position where it can be easily unplugged. The LED indicates the operational readiness of the power supply. The power indicator (11) is automatically activated when the device is switched on.

To put the power supply out of operation just disconnect the power supply from the mains supply.

 Do not pull the cable to disconnect the power adapter from the mains.

### Charging of the OMEGA 600

 Charge the device outside of the patient environment (at least 1.5 metres from the patient or patient support pursuant to IEC 60601-1, see Fig. 7).


Fig. 7

## Charging via USB

For charging, connect the USB power supply unit (e.g. E4-USB) to the USB-C socket (9) of the device.

 In charging mode, the LED lighting is switched off and the illumination of the device is no longer supplied with power.

 Do not let the power cable become taut as this could damage the device or pose a tripping hazard.

A description of how to start up and operate the E4-USB power supply unit is provided in a separate instruction for use.

## Charging via Wall Charger CW1

Use the Wall Charger CW1 to charge the battery CB1 (10) of the device. The blue LED on the Wall Charger CW1 indicates contact with the device. The charge status indicator (7) is located on the battery compartment (8) of the device.

 Make sure that no foreign objects come between the device and the Wall Charger CW1. When using the Wall Charger CW1 to charge the device, make sure that Wall Charger CW1 is not contaminated in any form.

The setting up and operation of the Wall Charger CW1 is described in a separate instruction of use.

## Charging via Charging Case CC1

Use the Charging Case CC1 to charge the battery CB1 (10) of the device. Turn the OMEGA 600 off. Remove the battery CB1 (10) from the battery compartment (8) and insert it into the Charging Case CC1. The charge status indicator of the Charging Case CC1 is located on the upper side of it.

 Make sure that no foreign objects come between the battery CB1 (10) and the Charging Case CC1. When using the Charging Case CC1 to charge the battery CB1 (10), make sure that the Charging Case CC1 is not contaminated in any form.

The setting up and operation of Charging Case CC1 is described in a separate instruction of use.

## Charge status indicator (7)

The charge status indicator (7) of the device is located on the battery compartment (8) of the device.

The charge status indicator (7) is automatically activated when the device is switched on.

Orange/green/green:	66–100 %
Orange/green:	33–66 %
Orange:	10–33 %
Orange flashing:	<10 %

## Operation

### Initial optical set up

Remove the protective dust cover (16) and place aside for reattaching after the examination. Unlock the adjustment lever (5) so that the optics unit (13) is free to move. Place the instrument on your head and adjust the height and circumference by means of the width adjustment (1) and height adjustment (2), respectively, until a comfortable fit is achieved. The rear part of the headband can be adjusted according to personal preference. Adjust the optics unit (13) to a position as close as possible to your eyes and centred to your face, then lock it in this position using the adjustment lever (5). Adjust the eyepieces (19) horizontally to match your own personal pupillary distance. The initially mounted eyepieces (19) incorporate +2D lenses that can be exchanged with neutral lenses (0D). Switch on the light by turning the brightness control (6) clockwise. You should now see the illumination spot being centred to your view at a distance of about 40 cm and you should be able to observe a pencil-sized object sharply focused at this distance. If you cannot focus on the pencil-sized object, you may have to adjust your distance to the object. Alternatively, you may try exchanging the initially mounted eyepieces (+2D lenses) (19) with the eyepieces containing neutral lenses (0D). More details on how to interchange the eyepieces (19) can be found in the section "Maintenance". If the illumination spot is not centred, you can adjust it vertically by twisting the illumination height adjustment (15) and you can adjust it horizontally by simply rotating the whole instrument slightly to the desired side. If a proper alignment has not been achieved, repeat steps as above. Correct adjustment of eyepieces (19) which match your own personal pupillary distance is particularly important when examining through small pupils. Each user should adjust the setting to match their own personal pupillary distance.

## Setting the brightness

Adjust the brightness by using the brightness control (6). To increase the brightness, turn the brightness control (6) clockwise. Once you reach a latch, you are operating in the highest brightness range in standard mode. If you turn the brightness control (6) further clockwise and over the latching mechanism, you switch into the visionBOOST setting. When the end stop is reached, the highest brightness range in the visionBOOST is operating.

To decrease the brightness or to turn the device off, turn the brightness control (6) counterclockwise. When the end stop is reached, the device is turned off.

It is recommended to start with the lowest possible brightness and then increase as required for the examination. The visionBOOST might be used for the examination of patients with media opacities, e.g. cataract.

## Aperture selection lever (14)

Three different aperture sizes and a diffuser aperture can be selected by means of the aperture selection lever (14). The choice of aperture depends mainly on the size of the patient's pupil. The diffuser is useful for the examination of the periphery.

## Filter selection lever (18)

In addition to the unfiltered option, the filter selection lever (18) can be used to select an interference red-free filter, a blue filter or a yellow filter. These are switched into the illumination beam.

The interference red-free filter can be used to view changes of the retina (e.g. new vessels or retinal nerve fiber layer defects) and highlights whitish portions of the retina, if present.

The blue filter can be used to do fluorescein angiography.

The yellow filter can be used to reduce discomfort for the patient and photochemical hazard by reducing blue light.

## Stereoscopic adjustment lever (17)

This feature ensures the best possible stereopsis for any pupil size and from any observation angle possible (e.g. peripheral view). When examining a dilated pupil head on, it is recommended to move the stereoscopic adjustment lever (17) into the forward position. When examining the periphery of the eye or in case of undilated pupils, it is recommended to move the stereoscopic adjustment lever (17) towards the side position (turn to your left). The stereoscopic adjustment lever (17) can be positioned anywhere between these two settings to select the optimum 3 dimensional view in any possible situation. The illumination beam can be adjusted vertically by the illumination height adjustment (15).

## Use of the flip-up function

While wearing the instrument, the optics unit (13) can be folded up into the rest position. To bring the optics unit (13) back into the working position, simply fold it down again. There is no need to readjust the instrument.

## Teaching mirror (TM)

The available teaching mirror can be attached to the optics unit (13) instead of the dust cover (16).

To be able to follow the examination of the main observer, the secondary observer stands next to the main observer and looks into the TM attached to the device from the side. The observation beam is split by the divider mirror, which approximately halves the brightness of the image for both the main and the secondary observer.

 The TM is intended for training purposes only.

## Hygienic reprocessing

The instruction is available:  
- online at [www.heine.com](http://www.heine.com)  
- in a paper version which you can request from the address listed

## Maintenance

There is no regular maintenance required. The following maintenance should be done outside of the patient environment (at least 1.5 metres from the patient or patient support pursuant to IEC 60601-1, see Fig. 7) and when deemed necessary.

### Changing the position of the brightness control (6)

The position of the brightness control (6) can be fitted on the right or left side onto the preferred socket for brightness control (4). To remove the brightness control (6), turn the device off and pull the dial out to gain access to the retaining screw. Remove the screw to detach the holder by using the enclosed offset screwdriver (see Fig. 8). Gently pull out the holder off the socket for brightness control (4) paying attention to the retaining clips. Remove the cover from the other socket for brightness control (4) and insert the holder by aligning the retaining clips. Insert the screw and tighten hand tight with the enclosed offset screwdriver (see Fig. 8). Reinsert the brightness control (6). Cover the spare socket for brightness control (4) with the cover. To activate the brightness control (6), remove the battery CB1 (10) from the OMEGA 600 or unplug the OMEGA 600 wired. Once power is reapplied the brightness control (6) will be fully functional.

Fig. 8

### Change of paddings (3)

To change the paddings (3), carefully pull on the corresponding padding (3), which is attached to the device via a Velcro fastening. To attach the paddings (3) to the device, press the corresponding padding (3) onto the Velcro fastener.

### Change of eyepieces (19)

To remove the eyepieces (19) unscrew them counterclockwise until they detach. To attach the eyepieces (19) screw them clockwise until they are hand tight.

## Service

### Exchange of the battery CB1 (10)

⚠ Exchange the battery CB1 (10) outside of the patient environment (at least 1.5 metres from the patient or patient support pursuant to IEC 60601-1, see Fig. 7).

Turn the OMEGA 600 off and disconnect the USB power supply unit. Remove the battery CB1 (10) from the battery compartment (8) by gently pulling and pushing the mechanical interlock at the same time. Do not touch the contacts on the battery CB1 (10). After inserting a new battery (10) make sure that the mechanical interlock is engaged.

Removing and reinserting of the battery CB1 (10) will reset the system.

The battery CB1 (10) only needs to be replaced if it can no longer be charged to a sufficient capacity. This is usually the case when the operating time of the battery CB1 (10) is reduced.

For the rechargeable batteries (10) two-year guarantee, please tell us the serial number of both the newly installed rechargeable battery (10) and your OMEGA 600:

[www.heine.com/OMEGA600/battery-change](http://www.heine.com/OMEGA600/battery-change)

## General Notes

⚠ The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

The expected life cycle amounts with designated use and the observation of warning and safety information as well as the maintenance instructions up to 7 years. Beyond this period, the product may continue to be used if it is in a safe and good condition.

Note to the user and/or patient:

All serious incidents that occur in connection with the product must be reported to HEINE Optotechnik GmbH & Co. KG and the member state's competent authority.

Charge your device during longer storage periods to protect the battery from discharge or store the battery and the device separate from each other.

Let the device accommodate to ambient conditions when it is cold.

Change the ocular lenses only in a clean environment in order to avoid dust to enter the indirect ophthalmoscope.

During non-use or storage of the device, cover the observation optics with the dust cover (16) to avoid dust on the exterior side of the glass.

Only use the indirect ophthalmoscope when the lenses are clean.

## General Warnings

⚠ Use only CE marked USB power supplies (5 V) from reputable manufacturers that conform to the safety requirements from IEC 60601-1 Medical electrical equipment.

Do not use power supplies where damage has been detected.

Check the correct operation of the device before use! Do not use the device if there are visible signs of damage or the light begins to flash.

Do not use the device in presence of flammable gases / liquids, or in an oxygen rich environment.

This product is not allowed to enter or be used in areas with strong magnetic fields e.g. MRI scanners.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

Do not use the device including its accessories and options outdoors.

Heating during operation is normal and harmless.

Make sure that the device is not exposed to direct sunlight.

For examination of the posterior segments of the eyes, use the indirect ophthalmoscope only in combination with ophthalmoscopy lenses.

## Light exposure hazard

 Because prolonged intense light exposure can damage the retina, the use of the device for ocular examination should not be unnecessarily prolonged, and the brightness setting should not exceed what is needed to provide clear visualization of the target structures. The retinal exposure dose for a photochemical hazard is a product of the radiance and the exposure time. If the value of radiance were reduced in half, twice the time would be needed to reach the maximum exposure limit.

While no acute optical radiation hazards have been identified for direct or indirect ophthalmoscopes, it is recommended that the intensity of light directed into the patient's eye be limited to the minimum level which is necessary for diagnosis. Infants, aphakes and persons with diseased eyes will be at greater risk. The risk may also be increased if the person being examined has had any exposure with the same instrument or any other ophthalmic instrument using a visible light source during the previous 24 hours. This will apply particularly if the eye has been exposed to retinal photography.

ISO 15004-2: Group 2 (LED)

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater the risk of ocular damage. An exposure time with this instrument at maximum intensity of longer than (see table exposure guideline) will lead to the guideline hazard value being exceeded.

ANSI Z80.36-2016: Group 2 (LED)

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater is the risk of ocular damage. Exposure to light from this instrument when operated at maximum intensity will exceed the recommended maximum exposure (RME) of 2.2 J/cm<sup>2</sup>, unless additional action is taken by the user to minimize exposure, after \_\_\_ min (see table exposure guideline). The risk of retinal injury at an exposure of 2.2 J/cm<sup>2</sup> is not high, but because some patients may be more susceptible than others, caution is advised if this radiant exposure value is exceeded. However, because of a significant risk of injury at exposures exceeding 10 J/cm<sup>2</sup>, the user should avoid exposures longer than \_\_\_ min. (see table exposure guideline).

## Exposure Guidelines

Data for the highest brightness setting in visionBOOST

	Working distance	Duration according to ANSI Z80.36-2016	Duration according to EN ISO 15004-2:2007
Without ophthalmoscopes	400 mm *	30 sec	2 min 44 sec
With ophthalmoscopes***	400 mm **	28 min	169 min

\*) from instrument

\*\*) from instrument to HEINE A.R. Aspheric Ophthalmoscopy Lens (A.R. 16D), diameter: 54 mm, focal length: 16 dpt.

\*\*\*) HEINE A.R. Aspheric Ophthalmoscopy Lens (A.R. 16D), diameter: 54 mm, focal length: 16 dpt.

## Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

 Dispose the rechargeable battery (10) at your local collection point.

The appendix contains following tables

- Electromagnetic disturbances – Requirements and tests
- Technical specification
- Explanation of the used symbols

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## Electromagnetic disturbances – Requirements and tests

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environments.

Statement for the operational environments	<p>Inside professional healthcare facilities except for: near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances are high.</p> <p>The supply voltage quality should be that of a typical hospital environment.</p> <p>Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.</p>
Performance features of the ME system that have been determined to be essential to the performance	None
Necessary instructions for maintaining basic safety and essential performance with regards to electromagnetic disturbances for the expected life cycle	
Warning	<p>Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.</p> <p>Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.</p> <p>Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</p>
Note	<p>Flickering of the LED of the OMEGA 600 wired is possible because of radiated electromagnetic fields.</p> <p>Brief illumination of the LED during charging of the OMEGA 600 is possible because of electrostatic discharge.</p>
A list of all cables, transducers and other accessories that are relevant for the EMC compliance	EMC compatibility is only ensured if original HEINE spare parts, accessories and power sources are used as described in the chapter „Accessories“. The EMC compatibility when using power sources from other manufacturers must be evaluated by the user.
<b>Test</b>	<b>Compliance</b>
RF emissions CISPR 11	Group 1 Class B
Conducted emissions (EN 55011/CISPR 11)*	Passed
Radiated emissions (EN 55011/CISPR 11)	
Harmonic current emissions (IEC 61000-3-2)*	
Voltage changes, voltage fluctuations and flicker (IEC 61000-3-3)*	
Immunity	See attached immunity test levels

## Immunity test levels

Test	Test level	
	IEC 60601-1-2 test levels	Compliance test levels
Electrostatic Discharge (IEC 61000-4-2)	Contact Discharge: $\pm 8$ kV Air Discharge: $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV	
Radiated RF EM fields (IEC 61000-4-3)	3 V/m 80–2700 MHz 80 % AM at 1 kHz	
Electrical fast transients   bursts (IEC 61000-4-4)*	$\pm 2$ kV 100 kHz repetition frequency	
Surges (IEC 61000-4-5)*	$\pm 0.5$ kV, $\pm 1$ kV	
Conducted disturbances induced by RF fields (IEC 61000-4-6)*	3 V 0.15 MHz – 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz, FM: $\pm 5$ Hz deviation: 1 kHz sine; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m	
Power frequency magnetic fields (IEC 61000-4-8)	30 A/m; 50Hz or 60 Hz	30 A/m; 60 Hz
Voltage dips (IEC 61000-4-11)*	0 % $U_T$ ; 0.5 cycle; at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % $U_T$ ; 1 cycle and 70 % $U_T$ ; 25/30 cycles Single phase: at 0°	
Short interruptions (IEC 61000-4-11)*	0 % $U_T$ ; 250/300 cycles	

\*n/a: "Not applicable" in the internally powered mode

**Technical Specification OMEGA 600 and OMEGA 600 wired**

Environmental conditions for operation	+10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Classification according to EN ISO 15004-2:2007 and ANSI Z80.36-2016	Group II The classification was performed together with a Ø 54mm 16 Diopter HEINE ophthalmoscopy lens.

**Technical Specification OMEGA 600**

CB1	Li-Po cell
Input	USB 2.0 Type C: 5 V, 1.2 A
Power consumption	6 W
Protection class	Charging: class II Operating: internally powered
Charging time	typ. 1.5 h
Operating time (at maximum charge capacity)	typ. 4 h
Operating time (visionBOOST)	typ. 1.5 h
Weight	475 g incl. rechargeable battery

**Technical Specification OMEGA 600 wired**

Input	100–240 V~   50–60 Hz   160–80 mA
Power consumption	6 W
Protection class	Class II
Weight	655 g

**Accessories**

CW1 – Wall Charger	X-095.17.320
CC1 – Charging Case	X-000.99.091
E4-USBC (USB-C cord approx. 2 m)	X-000.99.300

**Options**

TM2 – Teaching Mirror	C-000.33.212
OMEGA 600 Breath Shield	C-000.33.019
Fundus charts Pad with 50 pcs.	C-000.33.208

**Spare parts**

CB1 – OMEGA 600 Battery	X-007.99.687
OMEGA 600 eyepiece +2D	C-000.17.116
OMEGA 600 eyepiece 0D	C-000.17.115

**Erläuterung der verwendeten Symbole**  
**Explanation of utilized symbols**  
**Explication des symboles utilisés**  
**Explicación de los símbolos utilizados**  
**Spiegazione dei simboli utilizzati**  
**Förklaring av symboler som används**  
**Verklaring van de gebruikte symbolen**  
**Forklaring af de anvendte symboler**  
**Symbolforklaring**  
**Käytettyjen symbolien selitys**  
**Explicação dos símbolos utilizados**

<b>CE</b>	<p>CE-Kennzeichnung kennzeichnet die Übereinstimmung mit der Europäischen Medizinprodukteverordnung (EU) 2017/745.</p> <p>The CE mark indicates that the product complies with the medical device regulation (EU) 2017/745.</p> <p>Le marquage CE indique que le produit est conforme ou au règlement (UE) 2017/745 relatif aux dispositifs médicaux.</p> <p>La marca CE indica que el producto cumple el Reglamento europeo sobre productos sanitarios (UE) 2017/745.</p> <p>Il marchio CE indica la conformità con il regolamento sui dispositivi medici (UE) 2017/745.</p> <p>CE-märkning markerar en överensstämmelse med förordningen om medicinska produkter (EU) 2017/745.</p> <p>CE-markering duidt de overeenstemming aan met de verordening betreffende medische hulpmiddelen (EU) 2017/745.</p> <p>CE-mærkningen angiver overensstemmelse med forordningen om medicinsk udstyr (EU) 2017/745.</p> <p>CE-merket angir at produktet er i samsvar med forskriften om medisinsk utstyr (EU) 2017/745.</p> <p>CE-merkintä tarkoittaa, että laite lääkinnällisiä laitteita koskevan asetuksen (EU) 2017/745 kanssa.</p> <p>O símbolo CE identifica a concordância com o regulamento de Dispositivos Médicos (EU) 2017/74.</p>
	<p>Katalog- oder Bestellnummer          Catalogue- or order number          Numéro de catalogue ou de commande          Número de catálogo o de pedido          Codice catalogo e di dell'ordine numero          Katalog- eller Beställningsnummer          Catalogus- of bestelnummer          Katalog- eller Ordrenummer          Katalog- eller bestillingsnummer          Luettelo- tai viitenumero          Número de catálogo ou pedido</p>

<b>REF</b>	<p>Katalog- oder Bestellnummer          Catalogue- or order number          Numéro de catalogue ou de commande          Número de catálogo o de pedido          Codice catalogo e di dell'ordine numero          Katalog- eller Beställningsnummer          Catalogus- of bestelnummer          Katalog- eller Ordrenummer          Katalog- eller bestillingsnummer          Luettelo- tai viitenumero          Número de catálogo ou pedido</p>
------------	---

	<p>Hersteller          Manufacturer          Fabricant          Fabricante          Produttore          Tillverkare          Fabrikant          Producent          Produsent          Valmistaja          Fabricante</p>
---	--

	<p>Hersteldatum          Date of manufacture          Date de fabrication          Fecha de fabricación          Data di produzione          Tillverkningsdatum          Productiedatum          Produktionsdato          Produksjonsdato          Valmistuspäivä          Data de fabricação</p>
---	---

	<p>Getrennte Sammlung von Elektro- und Elektronikgeräten. (Europäische WEEE Richtlinie)</p> <p>Product bearing this symbol may not be disposed of together with general household waste, but instead requires separate disposal according to local provisions. (European Waste Electrical and Electronic Equipment Directive, WEEE)</p> <p>Tri sélectif des appareils électriques et électroniques. (Directive européenne DEEE)</p> <p>Desechado separado de aparatos eléctricos y electrónicos. (Directiva Europea RAEE)</p> <p>Raccolta differenziata di apparecchi elettrici ed elettronici (direttiva europea RAEE).</p> <p>Separat insamling av elektriska och elektroniska apparater (det europeiska WEEE-direktivet).</p> <p>Gescheiden inzameling van elektrische en elektronische apparaten (Europese AEEA-richtlijn).</p> <p>Separat innsamling af elektrisk og elektronisk udstyr (det europæiske WEEE-direktiv).</p> <p>Produktet med dette symbolet skal ikke kasseres sammen med vanlig husholdningsavfall, men krever separat kassering i henhold til lokale bestemmelser. (European Waste Electrical and Electronic Equipment Directive, WEEE)</p> <p>Sähkö- ja elektroniikkalaitteille tarkoitettu erillinen keräyspiste (eurooppalainen WEEE-standardi).</p> <p>Coleção separada de aparelhos elétricos e eletrônicos (Diretrizes Europeias WEEE).</p>
---	--

	<p>Batterien müssen einer zentralen Sammelstelle zugeführt werden. Dispose the rechargeable batteries at your local collection point. Mise au rebut au point de collecte local. Desechar en un punto limpio de la comunidad. Smaltimento in un centro di raccolta comunale. Avfallshanteras på ett kommunalt insamlingsställe. Afvoer naar een gemeentelijke inzamelocatie. Bortskaffelse på et kommunalt indsamlingssted. Avhænd hos ditt lokale afvalshåndteringsanlegg. Laite on hävitettävä paikalliseen keräyspisteeseen. Proceda à eliminação do aparelho em um ponto de coleta municipal.</p>
---	--

	<p>Zulässiger Temperaturbereich in °C für Lagerung und Transport          Temperature limits in °C for storage and transport          Plage de température admise en °C pour le stockage et le transport          Rango de temperatura permitida en °C para almacenar y transportar el producto          Temperatura ammessa in °C per conservazione e trasporto          Tillåtet temperaturintervall i °C för lagring och transport          Toegestane temperatueren in °C voor opslag en transport          Tilladt temperaturområde i °C ved opbevaring og transport          Temperaturbegrensning i °C for oppbevaring og transport          Näyttää pakkauksen sallitun säilytys- ja kuljetuslämpötilan (°C).          Limite de Temperatura permitida em °C para armazenamento e transporte</p>
---	--

	<p>Zulässiger Temperaturbereich in °F für Lagerung und Transport          Temperature limits in °F for storage and transport          Plage de température admise en °F pour le stockage et le transport          Rango de temperatura permitida en °F para almacenar y transportar el producto          Temperatura ammessa in °F per conservazione e trasporto          Tillåtet temperaturintervall i °F för lagring och transport          Toegestane temperatueren in °F voor opslag en transport          Tilladt temperaturområde i °F ved opbevaring og transport          Temperaturbegrensning i °F for oppbevaring og transport          Näyttää pakkauksen sallitun säilytys- ja kuljetuslämpötilan (°F)          Limite de Temperatura permitida em °F para armazenamento e transporte</p>
---	---

	<p>Zulässige Luftfeuchtigkeit für Lagerung und Transport          Humidity limitation for storage and transport          Humidité admise pour le stockage et le transport          Humedad del aire permitida para almacenar y transportar el producto          Umidità atmosferica ammessa durante il trasporto e la conservazione          Tillåten luftfuktighet för transport och lagring          Toegestane luchtvochtigheid voor opslag en transport          Tilladt luftfugtighed ved opbevaring og transport          Fuktighetsbegrensning for oppbevaring og transport          Sallittu ilmakestheus kuljetuksen ja varastoinnin aikana          Umidade do ar admissível para o armazenamento e transporte</p>
---	--

	<p>Zulässiger Luftdruck für Lagerung und Transport          Pressure limitation for storage and transport          Pression atmosphérique admise pendant le transport et le stockage          Presión de aire permitida para almacenar y transportar el producto          Pressione atmosferica ammessa durante il trasporto e la conservazione          Tillåten lufttryck för lagring och transport          Toegestane luchtdruk voor opslag en transport          Tilladt lufttryk ved oppbevaring og transport          Trykkbegrensning for oppbevaring og transport          Sallittu ilmanpaine kuljetuksen ja varastoinnin aikana          Pressão do ar admissível para o armazenamento e transporte</p>	<p>Nur in geschlossenen Räumen benutzen.          For indoor use only.          Utiliser uniquement dans des locaux fermés.          Sólo utilizar en espacios cerrados.          Utilizzare solo in ambienti chiusi.          Får endast användas i slutna rum.          Får sluitend in afgesloten ruimten toepassen.          Må udelukkende benyttes i lukkede rum.          Skal kun brukes i lukkede rom.          Käyttöön on sallittu ainoastaan suljettuissa tiloissa.          Utilizar apenas em espaços fechados.</p>
	<p>Vorsicht Bruchgefahr!          Fragile, handle with care!          Fragile ! Manipuler avec soin          Atención. Frágil.          Attenzione: pericolo di rottura!          Försiktigt! Risk för brott          Voorzichtig, kans op breuk!          Forsigtig, risiko for brud!          Ømtålig, behandles forsiktig!          Varo särkymisvaara!          Perigo de quebra!</p>	<p>          The Regulatory Compliance Mark (RCM)          The Regulatory Compliance Mark (RCM)          Marque réglementaire de conformité (RCM)          The Regulatory Compliance Mark (RCM)          The Regulatory Compliance Mark (RCM)          The Regulatory Compliance Mark (RCM)          The Regulatory Compliance Mark (RCM)          The Regulatory Compliance Mark (RCM)          Vaatimustenmukaisuusmerkki (RCM)          (RCM)- Marca de Conformidade Regulamentar</p>
	<p>Trocken lagern!          Keep dry!          Conserver au sec !          Conservar en un lugar seco!          Evitare ambienti umidi!          Förvaras torr!          Droog bewaren!          Opbevares tørt!          Hold tørt!          Säilytetään kuivassa paikassa!          Armazenar em ambiente seco!</p>	<p>          Geräte der Schutzklasse II          Class II equipment          Appareil de classe de protection II          Aparato de clase de protección II          Apparecchio di classe di protezione II          Apparat i skyddsklass II          Apparaat van beschermingsklasse II          Apparater i sikkerhedsklasse II          Klasse II-utstyr          Suojausluokan II laite          Aparelho da classe de proteção II</p>
	<p>Gebrauchsanweisung verbindlich befolgen.          (Hintergrundfarbe: blau, Vordergrundfarbe: weiß)          Follow instructions for use!          (Background color: blue, foreground color: white.)          Suivre le mode d'emploi.          (Couleur de fond : bleu ; couleur du premier plan : blanc)          Seguir obligatoriamente las instrucciones de uso.          (Color de fondo: azul, color de primer plano: blanco)          Attenersi obbligatoriamente alle istruzioni per l'uso.          (Colore dello sfondo: blu, colore in primo piano: bianco)          Bruksanvisningen ska alltid följas.          (Bakgrundsfärg: blå, förgrundsfärg: vit)          De gebruiksaanwijzing is bindend en dient gevolgd te worden.          (achtergrondkleur: blauw, voorgrondkleur: wit)          Følg alltid bruksanvisningen.          (Baggrundsfarve: Blå; forgrundsfarve: Hvid)          Følg bruksanvisningen!          (Bakgrunnsfarge: blå, forgrunnsfarge: hvit)          Käyttöohjetta on noudatettava tarkasti.          (Taufaväri: sininen, etualan väri: valkoinen)          Siga as instruções de uso!          (Cor de fundo: azul, cor de primeiro plano: branco)</p>	<p>          Wechselstrom (AC)          Alternating current (AC)          Courant alternatif (CA)          Corriente alterna (CA)          Corrente alternata (AC)          Växelström (AC)          Wechselstrom (AC)          Vekselstrøm (AC)          Vekselstrøm (AC)          Vaihinvirta (AC)          Corrente alternada (CA)</p>
	<p>Gleichstrom          Direct current DC          Tension continue          Tensión continua (CC)          Tensione continua          Likspänning          Gelijkspanning          Jævnspænding DC          Likestrøm dc          Tasajännite          Tensão continua (CC)</p>	<p>          Gleichstrom          Direct current DC          Tension continue          Tensión continua (CC)          Tensione continua          Likspänning          Gelijkspanning          Jævnspænding DC          Likestrøm dc          Tasajännite          Tensão continua (CC)</p>
	<p>Unique Device Identification</p>	<p>          Ladestandanzeige          Charge status indicator          Indicateur de l'état de charge          Indicador del estado de carga          Indicatore dello stato di carica          Indikator för laddningsstatus          Indicatielampje oplaadstatus          Ladestatusindikator          Ladestatusindikator          Lataustilan merkkipalo          Indicador do estado de carga</p>
	<p>Medical Device</p>	<p>5V ≈ 1.2A</p>
	<p>Gebrauchsanweisung          Instructions for use          Mode d'emploi          Manual de instrucciones          Istruzioni per l'uso          Bruksanvisning          Gebruiksaanwijzing          Bruksanvisning          Bruksanvisning          Bruksanvisning          Käyttöohjeet          Instruções de utilização</p>	

 **HEINE Optotechnik GmbH & Co. KG**  
Dornierstr. 6, 82205 Gilching, Germany  
[www.heine.com](http://www.heine.com)

V-200.00.328