

Hygienic reprocessing

HEINE mini-f

General warning and safety information



WARNING! This symbol indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent.



NOTE! This symbol is used for information regarding installation, operation, maintenance or repairs that are important but are not associated with dangers.

	Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines. The described reprocessing measures do not replace the specific rules applicable for your institution/ department.
	Hold the device in such a way that you do not accidentally shine the light into your eyes whilst using it. In the event of suspected contamination, carry out hygienic reprocessing of the instrument. Wear gloves when preparing and using the device. HEINE Optotechnik GmbH & Co. KG only approves the agents and procedures listed in this instruction. Hygienic reprocessing is to be carried out by persons with adequate hygienic expertise. Please consider the instructions of the manufacturer for the applied reprocessing media. Avoid liquid entry into the device.
Limitations on reprocessing	Periodically check the integrity of the device.

Wipe disinfection

1. Preparation

Switch off the device and let it cool down.

2. Wipe disinfection



Pay attention that all contaminated surfaces are completely moistened for the complete exposure time specified by the disinfectant manufacturer. If necessary, increase the number of wiping procedures and/or the number of wipes.



Preparation:

- When cleaning the device, hold it facing downwards to prevent liquids from entering it.
- Set the button① to the OFF position before cleaning.

Equipment:

- Disinfectant:
Super Sani-Cloth® by PDI (active ingredient based on alcohol and quaternary ammonium compounds)
- Nylon brush, soft bristles, length 15–40 mm
- Cotton buds

Implementation:

- Clean and disinfect the outer surfaces by hand using a damp cloth in circular motions for 1 minute.
- Pay particular attention to hard-to-reach areas.
- Dispose of the cloth.
- Take a new cloth and wipe again in circular motions for 1 minute.
- If there is heavy soiling around the button① it may be temporarily moved to the ON position.
- Whilst the surface is still damp, use a soft nylon brush to reach hard-to-reach areas such as the button① or ridges②.
- Whilst the surface is still damp, use a cotton bud for the clip area below the button①. To do this, gently lift the clip with your other hand.
- Take a new cloth and wipe again for 1 minute using circular movements.
- Dispose of the cloth.
- Remove any residue with a cloth moistened with deionised water and dry the device in accordance with the disinfectant manufacturer's instructions.



3. Inspection and function testing



- Check the device for any visible contaminants or abrasions. Reprocess again if necessary.
- Do not use the device if there are visible signs of damage..

4. Storage

Store it in such a way that it is protected from recontamination, dust and moisture.