



HEINE OMEGA 600 VET

indirect ophthalmoscope

A QUICK LOOK AT THE MOST IMPORTANT FEATURES.

We thank you for choosing HEINE OMEGA 600 VET and wish you a lot of success and fun working with your new instrument.



The video "How to set up your OMEGA 600" explains how to adjust the instrument to your personal needs.



HEINE OMEGA 600 VET

indirect ophthalmoscope



Flip-up function of the optics unit
For more convenient patient interactions



Eyepieces: To best accommodate your view, exchange the initially mounted eyepieces +2D with the neutral lenses 0D (supplied with the instrument)



Width adjustment



Teaching mirror TM2



Dust cover for front glass

If you want to share your experiences or discoveries, we would be highly pleased to see it:
#heine #iuseheine #OMEGA600



www.heine.com

