

# HEINE BETA X Retinoscope

Compact in size. Big on performance – easy to use, hard to beat.

The compact, lightweight HEINE BETA X Retinoscope features thoughtfully designed operating functions. These include precise and easy adjustment of the parallel beam path during examination, fine adjustment of the illumination level and a practical focusing ring for accurate reading of the cylinder axis.



# HEINE BETA X Retinoscope



## Focusing ring with integrated scale

- Fast, accurate cylinder axis reading
- Housing design enables quick and easy 90° rotation of the light beam via tactile feedback

## Optimised for hygienic reprocessing

- Scratch-resistant window made of tempered glass
- smoothSURFACE and recessed USB-C connection
- Dustproof optics with anti-reflective glass

## Durable, lightweight and ergonomic design

- Solid integrated aluminium frame
- Housing is thermo-moulded out of a special high-impact material blend optimised for hygienic reprocessing
- Precisely balanced for excellent ergonomics and secure handling
- Compact and lightweight design



**HEINE Optotechnik GmbH & Co. KG**

Dornierstr. 6, 82205 Gilching, Germany  
www.heine.com



We reserve the right to change  
specification without notice.

09|26. A-000.00.276 EN

**QUALITY  
MADE IN GERMANY**

