



HEINE iC1 Digital Dermatoscope.

Digital dermoscopy, redefined.

Sharp, sharper, sharpest. The HEINE iC1 is a device designed for capturing clinical and dermatoscopic images for documenting and supporting differential diagnosis of melanomas and other skin lesions. Developed in our unrivalled HEINE quality, the iC1 delivers sharp, shadow-free images. Thanks, to both our precise, aspherical optical system and, the 12-megapixel image quality when using an iPhone* 6s or later model. Download the exclusive HEINE Derm App, insert your iPhone and experience digital dermatoscopy like never before.

- ▶ **Toggle function for switching between polarised and non-polarised illumination mode**
- ▶ **Free HEINE DERM app**
- ▶ **LED in HEINE Quality – or LED^{HQ}**



iC1 digital dermatoscope

with contact plate with scale, mounting case smartphone iC1/5 for Apple iPhone 5s* and SE 1st generation*, USB cord with medical approved plug-in power supply	K-271.28.305
as above, but with mounting case smartphone iC1/6 for Apple iPhone 6/6s*	K-270.28.305
as above, but with mounting case smartphone iC1/7 for Apple iPhone 7*, 8* and SE 2nd generation*	K-273.28.305
as above, but with mounting case smartphone iC1/12 Pro for Apple iPhone 12 Pro*	K-277.28.305

* The Apple iPhone is not included in the scope of delivery. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.



- ▶ **Find all the details, accessories and order information at the QR code or at www.heine.com/iC1-product**

HEINE Optotechnik GmbH & Co. KG

Dornierstr. 6, 82205 Gilching, Germany
www.heine.com



We reserve the right to change specification without notice.

02|25. A-000.00.269 EN

**QUALITY
MADE IN GERMANY**

