

# HEINE DELTA X

## Dermatoscope

The entry point to professional standards. Everything you need for dermatoscopy – in a compact design.

### Intuitive and ergonomic

- Small instrument head and proven contact plate design – facilitates access to the examination area
- Ergonomic handling thanks to the angled handle
- Compact – ideal for mobile use

### Toggle function

Alternate between polarised and non-polarised light with the push of a button  
More information: [heine.com/en/Toggle](https://heine.com/en/Toggle)

### Additional light

For viewing the skin before the dermatoscopic examination

### Smart power source options

- Charge status indicator
- Automatic switch-off function
- Rechargeable with USB-C via cable or charging station

### Extremely sharp image, selectable polarisation eliminates reflections and glare

- Crystal clear image due to achromatic HEINE optical system
- Homogeneous, bright exposure with true colour rendering thanks to LED<sup>HQ</sup>
- 10x magnification
- Field of view 14mm (lens 16mm) Ø
- Focusing optics (dioptrre correction)

### Options for digital documentation

- Easy capture of digital images with the smartphone adapter (optionally available) for documentation or for a second opinion
- Convenient and intuitive image and data management with the HEINE DERManager



# HEINE DELTA X

## Dermatoscope

### Additional light for an overview before the dermatoscopic examination

- For viewing the skin before the dermatoscopic examination through the lens
- Homogeneous, bright illumination with true colour rendering

### smooth SURFACE

Safe, hygienic reprocessing thanks to smoothSURFACE and recessed connections



Product data

### Contact plates for multi or single use

- Contact plates made of high-quality optical polymer – for safety in hygienically problematic applications by replacing the contact plate
- Secure hold of the contact plate with the bayonet lock
- Quick change from contact to non-contact examination – simply remove the contact plate

### Optional: Contact plate with scale made of scratch-resistant glass

- For long-term use



### Durable and lightweight thanks to optimal material selection

- Solid integrated aluminium frame
- Housing is thermo-moulded out of a special high-impact material blend optimised for hygienic reprocessing
- Lightweight and precisely balanced
- LED<sup>HQ</sup>: Virtually unlimited durability and constant performance due to special thermal management (aluminium heat sinks)

HEINE Optotechnik GmbH & Co. KG

Dornierstr. 6, 82205 Gilching, Germany  
www.heine.com



We reserve the right to change specification without notice.

10|25. A-000.00.275 EN



QUALITY  
MADE IN GERMANY

